



TOROBlue® One HS KOD Kit

Cat No. : KAO-201 Vol: 50µL×40 Reactions

DESCRIPTION

TOROBlue® One HS KOD Kit includes 2× One HS KOD Mix with dye is a ready-to-use premixed solution containing genetically modified hotstart KOD DNA polymerase, PCR buffer, dNTPs, elongation accelerator and loading dye. This mix which enables fast PCR, which has an extension time of 5 sec/ kb by applying genetically modified KOD DNA polymerase and a new elongation accelerator. In addition, this mix provides greater efficiency and elongation capabilities than conventional PCR enzymes. In particular, these show greater amplification success from crude specimens.

FEATURES

- Fast:** can amplify the targets using the following very short conditions: 0~10 kb: 5 sec/ kb ; ≥10 kb: 10 sec/ kb
- Easy-to-use:** premixed with loading dye for directly loading the gel of Post-Staining or Precasting .
- High Fidelity:** exhibits approximately 80-fold higher fidelity than Taq DNA polymerase for cloning.
- Direct PCR:** effective for amplification from crude samples.

APPLICATIONS

- Direct PCR
- Colony PCR
- Amplification of NGS libraries
- Site direct gene mutation

COMPONENTS

The kit includes the following reagents, which can be used for 40 reactions for a total 50ul reaction volume, or 100 reactions for a total 20ul reaction volume.

Cat No. : KAO-201

2× One HS KOD Mix with Dye 1ml ×1 tube

PROTOCOL

1. Standard reaction setup

(1) Before preparing the mixture, all components should be completely thawed on ice.

Component	Reaction volume	Final Concentration
PCR grade water	XµL	
2× One HS KOD Mix with Dye	25µL	1 ×
10uM Forward primer	1.5µL	0.3µM
10uM Reverse primer	1.5µL	0.3µM
Template DNA	YµL	Genomic DNA: ≤200 ng / 50µL Plasmid DNA: ≤50 ng / 50µL cDNA: ≤750 ng / 50 µL Crude Sample: ≤5ul / µL
Total Volume	50µL	

(2) Gently mix the reaction solutions and spin down in microcentrifuge.

Notes:

- Primers should be 22-35 bases with $T_m \geq 65^\circ\text{C}$.
- Optimal primer concentration is 0.3µM. In the case of long targets (≥10 kb), reduced primers concentration (0.15µM) may give more effective amplification.
- When PCR yield is low, increased primers concentration (0.5 µM) may give more effective amplification.
- Decreased the amount of template DNA when no PCR Product.

2. Cycling conditions

[3-step cycle]

Predenaturation : 95°C, 30 sec.
 Denaturation : 98°C, 10 sec. ← 25-45 cycles
 Annealing : ($T_m - 5$) °C, 5 sec.
 Extension : 68°C, 5-10sec. /kb

[2-step cycle]

Predenaturation : 95°C, 30 sec.
 Denaturation : 98°C, 10 sec. ← 25-45 cycles
 Extension : 68°C, 5 -10sec. /kb

Notes:

- Extension time: 0-10kb, 5sec/kb ; ≥10 kb, 10sec/kb.
- Longer extension time may enhance efficiency. For amplification from a low copy DNA or crude sample, the extension time should be 10 sec. /kb.
- Poor amplification may be improved by changing the denaturation step to 94°C, 15sec.

3. Electrophoresis

- (1) Directly loading the PCR products and the suitable TOROBlue® Loading Marker into the Post-Staining or Precasting agarose gel.
- (2) Run on agarose electrophoresis to detect PCR products and marker.

RELATED PRODUCTS

→ Loading Marker:

BDL-029 TOROBlue® 100bp DNA loading Marker
 BDL-037 TOROBlue® 1.5kb DNA loading Marker
 BDL-036 TOROBlue® 5000bp DNA loading Marker
 BDL-031 TOROBlue® 1kb DNA loading Marker
 BDL-200 TOROBlue® DL2000 DNA loading Marker

→ Nucleic Acid Staining Dye:

RSD-105 TORORed® Nucleic Acid Staining Dye
 GSD-105 TOROGreen® DNA Staining Dye

→ Reagents for Electrophoresis :

AGR-100 Agarose
 TBE-054 5×TBE Buffer

→ Instruments for Electrophoresis :

MINI-300 Mini Pro 300V Power Supply
 MJ-105A Mini Horizontal Gel Electrophoresis System
 MBE-300 Mini BluView Transilluminator

STORAGE

This kit should be kept at -20°C.