

TOROGreen® 5G qPCR Premix 2.0

【Catalogue Number】 QST-200

【Packing Information】 1000 reactions for a total 20μL reaction volume.

【Description】

TOROGreen®5G qPCR Premix 2.0 is an 2×master mix for intercalator-based realtime PCR with SYBRGreen I, which contains all components except for the primer. TOROIVD®5G DNA polymerase, a mutant hot-start Tth DNA polymerase modified by antibodies, allows high specificity and sensitivity for high-speed PCR reactions. The improved polymerase and reaction mixture combination also enables a high resistance to PCR inhibitors. ROX is added into the premix and can be applied to the qPCR cyclers that require a passive reference dye. The premix is suitable for high-speed PCR and enables accurate detection and quantification of targets, making it possible to obtain highly reproducible and reliable realtime PCR results over a wide dynamic range.

【Feature】

-Rapid and highly sensitive

This premix can achieve the rapid and highly sensitive quantification of a low-copy targets with probes and be suitable for the quantification of DNA viruses or cDNA at a low level.

-Room-temperature stable:

The specially optimized PCR buffer make the mix very stable at room temperature. Therefore, the performance is not easily decrease during storing and shipping.

-Inhibitor tolerant

The unique proprietary formulation of this kit allows robust performance even in the presence of substances that can normally inhibit PCR, such as heparin, hematin, or EDTA.

-Wide dynamic range

The master mix demonstrates excellent reproducibility over a wide dynamic range and provides efficient amplification over 8 logs of sample input.

【Components】

QST-200 can be used for 1000 reactions for a total 20μL reaction volume.

Cat NO.	Components	Size
QST-200	TOROGreen® 5G qPCR Premix 2.0	1.25 mL ×8 tubes/ Kit

【Primer Design】

- Primer length: 18~30bp **- GC content of primer:** 40~60% **-Target length:** ≤ 200 bp (optimally, ≤ 150bp)

-Checking the primers:

-NTC tests can distinguish unintended amplification products of primer dimers from the intended PCR products in SYBR Green I reactions. NTC test is required to verify the each of primers for assessing the extent of primer dimers. The primer of NTC with C_q<40 should be redesigned.

-Prepare a dilution series with five or more dilutions of template DNA. Perform qPCR assay using the diluted DNA with the newly designed primers and draw a standard curve. Confirm that the PCR efficiency is between 95% and 105% and R² is equal to or greater than 0.99. If the PCR efficiency or R² are outside of these ranges, the primers concentration and reaction conditions should be optimized. If this does not improve the result, the primers should be redesigned.

【Template DNA】

-Genomic DNA: Purified DNA, which would be used for general PCR, is also suitable for real-time PCR. In the case of mammalian genomic DNA, 1~10 ng genomic DNA is sufficient for real-time PCR.

-cDNA: Reverse transcription reactions from total or poly (A)⁺ RNA may be used directly, or after dilution for realtime PCR. Before the reverse transcription reaction, it is essential to assess the extent of genomic DNA

contamination with no-reverse transcription control.If genomic DNA contamination affects the Cq values, it is essential to be eliminated by DNase treatment.

【Detection】

- This reagent can be used in general detection devices,not needing ROX such as: LineGene(bioer); LightCycler (Roche); iCycler iQ, CFX96(Biorad/MJ); Thermal Cycler Dice(Takara);
- This reagent with 1× ROX can also be used in detection equipment using passive reference, such as: ABI PRISM 7000, 7700, 7900 ,7300; Step One, Step one plus etc.(ABI) , ABI PRISM 7500, 7500Fast(ABI); Mx3000P, 3005P, MX4000,etc.(Agilent).

【Protocol】

1. Preparation of the reaction mix

- This premix should be fully thawed at room temperature in the bags, gently vortexed and briefly centrifuged.

Notes: Due to the high concentration stabilizer , there may be crystal precipitation in the premix , which can be used normally after being fully thawed at room temperature

- Purified DNA or RT reactions can be may be used directly or after dilution.
- In order to reducing the artificial error of sampling , design the plate layout and sampling method by the number of the templates and primer pairs. According to the following two situations , the total reaction is divided into two parts for premixing and loading in the a thin-walled qPCR tube or plate at room temperature.

Fore more genes and less samples in one plate

Components	20μL reaction ×n	Operation
TOROGreen®5G qPCR Premix 2.0	10μL×n	Premix and Loading
Template DNA Dilutions	2μL×n	
2μM Forward primer	4μL×n	Premix and Loading
2μM Reverse primer	4μL×n	

Fore more samples and less genes in one plate

Components	20μL reaction ×n	Operation
TOROGreen®5G qPCR Premix 2.0	10μL×n	Premix and Loading
8μM Reverse primer	1μL×n	
8μM Reverse primer	1μL×n	
Template DNA Dilutions	8μL×n	Premix and Loading

- Gently mix the reaction solutions and spin down in microcentrifuge.

2. Set up the cycling conditions

For Bio-Rad CFX96,etc.				
1	Pre-denaturation	95℃	1min	1cycle
2	Denaturation	95℃	3sec	40
	Annealing/ Extension	60℃	5 sec	cycles
Data collection should be performed at the extension step.				

For ABI StepOne Plus,etc.				
1	Pre-denaturation	95℃	1min	1cycle
2	Denaturation	95℃	3sec	40
	Annealing/ Extension	60℃	10sec	cycles
Data collection should be performed at the extension step.				

For Roche LightCycler 96 / 480 systems,etc.				
1	Pre-denaturation	95°C	1min	1cycle
2	Denaturation	95°C	10sec	40
	Annealing/ Extension	60°C	20sec	cycles
Data collection should be performed at the extension step.				

For ABI 7500/7300 etc.				
1	Pre-denaturation	95°C	1min	1cycle
2	Denaturation	95°C	10sec	40
	Annealing/ Extension	60°C	30sec	cycles
Data collection should be performed at the extension step.				

Notes:

- The annealing temperature can be set to 55~65°C, depending on the primer T_m value.
- The annealing time should be set for 5~20 seconds. Longer annealing time results increased efficiency, and a shorter time decreases non-specific amplification.
- Data collection step should be longer than 10 sec.

【Storage】

This reagent can be stored at 2-8°C for 24 months and protected from light.
For longer storage, this reagent should be kept at -20°C and protected from light.

【References】

Bustin SA, Benes V, Garson JA, etc,al. The MIQE guidelines: minimum information for publication of quantitative real-time PCR experiments.ClinChem.2009,Apr;55(4):611-22.

【Contact】



TOROIVD TECHNOLOGY COMPANY LIMITED.

Head Office: Building #C1, 880 JiangYang Nan Road, Baoshan,Shanghai,China. 200439.

Factory: Building #20, 888 Zhujiang Road, Rudong ,Jiangshu,China. 226400

Tel: +86-21-68030217

Mail: market@toroivd.com

http: //www.toroivd.com