

DESCRIPTION

TOROBlue® qRT Premix with gDNA Eraser enables fast cDNA synthesis for 2-step RT-qPCR and includes a rapid step that eliminates genomic DNA contamination without loss of input RNA. Genomic DNA is eliminated by treating the RNA templates with gDNA Eraser Mix for 2-5 minutes at 37°C. 5×qRT Premix contains primers and buffer optimized for highly efficient synthesis of short-chain cDNAs suitable for qPCR. The protocol is simple, and the reaction can be completed in 10-25 min.

COMPONENTS

The kit includes the following reagents, which can be used for 100 reactions for a total 10μl reaction volume.

Cat NO : RTQ-201

RNA Dilution Buffer	1ml/tube
gDNA Eraser Mix	200μl/tube
5×qRT Premix	200 μl/tube
RT Enzyme Mix	100μl /tube

Notes:

- RNA Dilution Buffer can be used to dilute and store RNA instead of RNase-free water.
- gDNA Eraser Mix contains Heat labile dsDNA Nuclease, RNase inhibitor, and reaction buffer.
- 5×qRT Premix contains oligo dT primer, random primer, and reaction buffer.
- RT Enzyme Mix contains the highly efficient reverse transcriptase and an RNase inhibitor.

PROTOCOL

1. Preparation of the reagent and templates

- This kit should be fully thawed before use. Gently vortexed and briefly centrifuged.
- Purified template RNA can be may be used directly or after dilution.

2. Genomic DNA elimination reaction

- Prepare the following reaction in a thin-walled PCR tube on the ice.

Component	Volume
gDNA Eraser Mix	2 μL
Total RNA	X μL
RNA Dilution Buffer	5-XμL
Total	7μL

- Incubate at 37°C for 5 min. Store the reaction tube on ice.

Notes:

- The indicated time of gDNA elimination reaction is between 2-5min at 37 °C.
- Up to 0.5 μg of total RNA template for qPCR assay.

3. Reverse-transcription

- Transfer 2 μL 5 × qRT Premix and 1 μL RT Enzyme Mix into the above reaction tube. Gently vortexed and briefly centrifuged.
- RT reaction as the following condition.

37°C	15min
50°C	5min
98°C	5 min
- Store the reaction tube on ice.

Notes:

- The cDNA products can be used directly or after dilution for real-time PCR. Up to 20% of the synthesized cDNA solution can be added to the PCR reaction solution.
- The mutant M-MLV reverse transcriptase excels at high reaction temperatures (up to 60 °C). This step may increase the efficiency of the reverse transcription.

APPLICATION DATA

1. Efficiency of gDNA elimination.

RT Reagent: *TOROBlue[®] qRT Premix with gDNA Eraser (RTQ-201)*

Template: *Human Genomic DNA 50ng/reaction*

Experiment groups:

	4 × gDNA Eraser Mix	5 × qRT Premix	RT Enzyme Mix
A	(+)	(+)	(-)
B	(-)	(+)	(-)

qPCR Reagent: *TOROGreen[®] qPCR Master Mix (Code No.QST-100).*

Template: *cDNA 2 μ L/20 μ Lreaction*


Target:*GAPDH (Non cross-intron)*

Forward primer: *AAAAGGGCCCTGACAACTCTT*

Reverse primer: *GTCTTACTCCTTGGAGGCCA*

Instrument : *CFX-96 from Biorad*

Results:



No signal for the “A groups” indicates that 50ng gDNA was completely removed by gDNA Eraser.

2.Comparison of RT performance

RT Reagent: *TOROIVD[®] qRT Master Mix(RTQ-100).*

TOROBlue[®] qRT Premix with gDNA Eraser (RTQ-201)

Template: Template RNA:*MS2RNA (Roche)*

Forward primer: *GCCTTAGCAGTGCCCTGTCT*


Reverse primer:*AACATGCTCGAGGGCCTTA*

qPCR Reagent: *TOROGreen[®] qPCR Master Mix (Code No.QST-100).*

Template: *cDNA 2 μ L/20 μ Lreaction*

Instrument : *CFX-96 from Biorad*

Results:



Despite the gDNAcontamination, the results of RTQ-201 correlate highly with those of the RTQ-100. Both reagents showed highly linear standard curves in a broad concentration range.

STORAGE

This reagent should be kept at -20°C for 1 years.