




\*For research use only. Not intended for human or animal clinical diagnostic use

Portable, low-cost and easy-to-use qPCR cycler

# TOROIVD® 5G Reaction Mix

Cat NO. TS-2500 Size: 2.5mL×1



## Product Features

✓ **Easy-to-use:** dropping operation as Ag-test.

✓ **High compatibility:** compatible with all test tube of TOROIVD.

✓ **No contamination:** contains dUTP/UNG in the reaction mix.

✓ **High stability:** stored at room temperature for 1 months.

### TOROIVD TECHNOLOGY COMPANY LIMITED.

Head Office: Building #C1, 880 Jiang Yang Nan Road, Baoshan, Shanghai, China. 200439

Service Hotline: 400-688-2055

Website: [www.toroivd.com](http://www.toroivd.com)

Manufacturing Site: No.888, Zhujiang Road, Juegang Sub-district, Rudong County, Jiangsu Province

Email: [market@toroivd.com](mailto:market@toroivd.com)

